





FACE SHIELD AND GOGGLE

UV radiation is dangerous for unprotected eyes and skin. Users must protect themselves against UV radiation by wearing glasses or face shields. The MP-80 is recommended for the protection of the eye and the face. The MP-800 is a face shield with two lateral protections to cover the operator ears in addition to his eyes and face. Comfortable and efficient, the LP-70 glasses provide total protection for the eyes.

CONVERSION SCREEN

The conversion screen converts the 312 nm UV light into blue or white light. The FC-26.WL converts UV to white light and is ideal for autoradiographs or protein gels. The FC-26.Blue converts UV to blue light (480 nm) and it could be used for application such as GFP II, SYBR Safe®, SYBR Green® or SYPRO Orange®.



LAMP STANDS

The lamp stands or the handlers give you ease of use and versatility. Lamp handle and bracket let you conveniently mount your lamp below a horizontal surface. The lamp stand frees your hands.

SPECIFICATIONS

| Models | Type | Compatibility |
|--------|---------|-------------------------|
| SVL-30 | Stand | For VL-215 and VL-115 |
| SVL-6 | Stand | For VL-8, VL-6 and VL-4 |
| SMA | Handle | For VL-8, VL-6 and VL-4 |
| SMU | Bracket | For VL-8, VL-6 and VL-4 |